

For thorough extraction of stubborn upholstery soils.

Matrix **Impression** Upholstery Cleaning Concentrate

Upholstery soils tend to be much more stubborn than those found on carpeting. The presence of hair and body oil creates sticky binders that act like glue, holding dirt and grime to the fabric fibers.

Once these binders have been dissolved, it's essential that they are completely flushed away to prevent resoiling. For this cleaning challenge, the chemists at Matrix developed Impression Upholstery Cleaning Concentrate.

Impression is formulated to withstand extreme heat, so it works with hot water extraction to remove the soils, hair, and body oils that were broken down during pre-treatment and easily lift them away. Upholstery is left clean and residue free.





Dilution and Usage Instructions

pH: 11.0-13.0.

Dilution: 1-2 oz per gallon of water in portables and 1 gallon to 5 gallons of water for concentrated truck mount stock solution. Meter at 3 GPH.

Usage Instructions:

Upholstery cleaning is a more delicate cleaning procedure than general carpet cleaning. Upholstery cleaning chemicals are formulated to be hypoallergenic due to increased contact with skin. Products used on upholstery must always be tested. Do this by diluting the chemical at the dilution you would use to clean the entire fabric. Look for bleeding, crocking, and color transfer while testing. Check the threading and piping for any abnormalities or different colors. Open up the cushions at the zipper and look for any magic marker used to mark cushions. **Common symbols on tags for upholstery are X - vacuum only, S - solvent only, W - wet clean, and S/W - solvent or wet clean.**

1. Vacuum well. Be sure to use a crevice tool in the corners and crevices. Never stick your hands inside the cushion, as there is no telling what could be inside.
2. Using a pump-up sprayer, apply Outset Upholstery Pre-Spray.
3. Use a horsehair brush to agitate heavily soiled areas.
4. Allow 10 minutes of dwell time for soil suspension.
5. Using an extraction machine, extract furniture with Matrix Impression Upholstery Cleaning Concentrate.
6. Dry vacuum all excess moisture, especially around the arms and cushions.
7. Speed dry with air movers.

