

Rescue Mat System

(Model: F286)

Symbol Definition

- * - multiple components to choose from for different models: voltage, size, color, etc.
- + - if you order this item, you may need to order additional items marked with same number of +'s in order to upgrade or complete repair

	Part Description	Specifications/Model	Part #	Qty
3	Gasket	.375 x .125 Closed cell (sold in 50' Roll)	03-00201	117 ft
	Box	46 x 32 x 6 Rescue Mat	06-01464	1
	Foam	12 x 12 x 1/2 Eps 1 lb w/ 5 1/2 hole	06-01466	5
	Bag	bag - 12 X 12 poly 1.5 mil	06-01467	1
	Instructions	Instructions - Rescue Mat	07-00438	1
	Hose Cuff	1.5 Threaded x 2.0 Slip Vinyl Black	13-00063	6
	Sleeve	2 x 3 Vinyl Black	13-00076	2
	Hose Assembly	1.5 x 7' Vacuum	13-00090	3
	Fitting	2 x 2 x 2 Barbed Tee (Non Tapered)	13-01550	5

Rescue Mat Pictorial Drawing



NOTES:

1. APPLY GASKET (ITEM 3) EVEN WITH EDGE OF PANEL